






N-EVAP Replacement Parts List

Item Description	Part Number	Price	Picture
Connector tube for gas supply, 6ft coiled polyurethane plastic	NA1101	\$29.50	
Swivel tube fitting for stainless steel gas supply tube	P1204	\$16.20	
Paper filled air filter for N2 supply, Z-purge or air compressor	NA0403	\$10.90	
Thermometer, 6in long, 100°C for N-EVAP water bath	NA1110	\$38.00	
Thermometer, 6in long, 150°C for N-EVAP dry bath	NA1121	\$38.00	
Valve tube assembly chrome with rear mount hose barb, for N-EVAPs made after 1997	P0627	\$34.30	

Luer and washer in chrome plated brass, for N-EVAP valve tube	P0607	\$8.80	
Nylon 3-part tube fitting for N-EVAP valve tube	P0634	\$4.50	
Hoist pulley and band spring For 6 and 12 Position N-EVAPs For 24 Position and Auto-20 N-EVAPs For 34 and 45 Position N-EVAPs	P1224-M P1224 P1524	\$38.00 \$38.00 \$52.00	
Flow meter 0-10 LPM, for 6 and 12 Position N-EVAPs 0-20 LPM, for 24 Position N-EVAPs 0-30 LPM, for 34 and 45 Position N-EVAPs	NA1221 NA1521 NA1421	\$128.00 \$128.00 \$128.00	
Gas Supply Tube	P1230	\$28.00	
Thumb Screw	P0612-N	\$3.40	